
EN – USER MANUAL	2
CZ – UŽIVATELSKÝ MANUÁL.....	4
SK – NÁVOD NA POUŽITIE	7
PL – INSTRUKCJA OBSŁUGI.....	9
DE – BENUTZERHANDBUCH	11
FR – MANUEL D'UTILISATION	13
IT – MANUALE D'USO	15
NL –GEBRUIKERSHANDLEIDING	17
ES – MANUAL DE USUARIO.....	19
PT – MANUAL DO UTILIZADOR.....	21
HU – HASZNÁLATI ÚTMUTATÓ	23
HR – UPUTSTVO ZA UPORABU	25
SI – UPORABNIŠKI PRIROČNIK	27
RS – UPUTSTVO ZA UPOTREBU	29
RO – MANUAL DE UTILIZARE	31
BG – РЪКОВОДСТВО ЗА ПОЛЗВАНЕ	33
UA – ІНСТРУКЦІЯ КОРИСТУВАЧА.....	35
DK – BRUGERMANUAL	37
FI – KÄYTTÖOHJE.....	39
SE – ANVÄNDARMANUAL.....	41
GR – ΕΓΧΕΙΡΙΔΙΟ ΧΡΗΣΤΗ.....	43
LT – NAUDOJIMO INSTRUKCIJA.....	45
LV – LIETOŠANAS ROKASGRĀMATA.....	47
EE – KASUTAJAKÄSITLUS	49
TR – KULLANIM KILAVUZU.....	50



A

B

C

Fitdays App

EN – USER MANUAL

Package Contents: Smart scale + measuring handle, user manual

Technical Specifications: dimensions: 333 x 305 x 33 mm; maximum capacity: 180 kg (max. 400 lb); power supply: 4x AA batteries (not included); measurement method: BIA (Bioelectrical Impedance Analysis) using 8 electrodes; app: Fitdays

Device Description:

1. Thumb electrode on the handle, 2. Measuring handle, 3. LED display, 4. Tempered glass surface, 5. Stainless steel plate, 6. Handle electrode for four fingers, 7. Scale feet, 8. Bottom plastic base, 9. Battery compartment cover (reset hole inside)

Getting Started:

- **Battery Installation:** Insert 4xAA batteries into the scale and place the scale on a hard and level surface.
- **Downloading and Installing the App:** Scan the QR code or search for the "Fitdays" app in the App Store / Google Play. Download and install the app.
- **Registration and Account Setup:** Open the Fitdays app and register. Fill in personal details and confirm.

Note: Incorrectly entered data (gender, age, height, etc.) can affect the accuracy of results. **Athlete Mode:** Active athletes (more than 10 hours of intense training per week) should turn on "Athlete Mode" in their profile for a more accurate calculation of muscle mass.

- **Turning on Bluetooth:** Turn on Bluetooth on your phone.
- **iOS:** Allow the app access to Bluetooth.
- **Android:** Turn on Bluetooth, GPS and allow access to location.
- **Activating the Scale:** Lightly tap the scale with your foot to activate it. Wait until "0.0 kg" is displayed.
- **Pairing the Device:** Go to the app: Account -> Device -> Add Device + -> select **WG BodyExpert**.
- **Unit Switching and Reset:** Gently insert a thin object (e.g., a paperclip) into the reset hole (inside the battery compartment) and press briefly to switch units (kg/lb). Long-press for about 6 seconds; when "CLr" appears on the display, the product has been reset and you can restart the measurement.

Measurement:

- Place the scale on a hard and level surface (e.g., tiles, laminate flooring, Figure A).
- Activate the scale by lightly tapping it with your foot. Wait until the display shows "0.0" (kg/lb).
- Step on the scale barefoot – feet must be in full contact with the electrodes.
- Grasp the measuring handle with both hands.
- Place your thumbs on the front electrodes on the top side of the handle.
- Place the tips of your other four fingers on the back electrodes of the handle.
- Extend your arms in front of you pointing downwards so they form an angle of approximately 45° with your body and do not touch your hips (Figure B).
- Stand upright during the measurement with your legs slightly apart and ensure your thighs do not touch. If necessary, separate them with thin clothing or a sheet of paper to avoid distorting the results.
- Stay still during the measurement and do not move your arms or body. The measurement progress is visible in the app.
- To obtain accurate results, remain standing on the scale and wait until the measurement is completed.
- **Important:** When measuring in socks or shoes, the scale will only measure weight.
- **User Selection:** If the scale is used by multiple people, select your profile in the app before measuring. If the scale detects a significant difference from the previous record, it will alert you to check the account to avoid data confusion. Without a selected profile, the device may only display total weight.
- **Baby Mode:** Allows you to weigh a baby, child or small animal. **Activation:** In the app by tapping "+" in the top right corner and selecting "Baby Mode". **Procedure:** The parent first stands on the scale alone, then picks up the child – the scale automatically calculates the child's weight as the difference.

Offline Weighing:

The first connection with the app is necessary for the scale to function correctly. This activates your profile and body fat measurement; you do not need to have the app running for subsequent weighings.

- **2 kg Rule:** Offline measurements are saved to your profile only if your current weight does not differ from the last one by more than 2 kg. This is protection against user confusion.
- **Automatic Synchronization:** All data measured without a phone will be automatically added to the app during the next connection.

Special Feature: Voice Assistant

The scale supports a voice assistant function that guides you through the measurement process with voice instructions. The feature is available in the following languages: English, Polish, Italian, French, Chinese, and Vietnamese. Settings in the app: Account → Device → Voice Assistant (on/off) Below, you can select the assistant language.

The scale measures the following data:

- **Basic Data:** Weight, BMI (Body Mass Index), Ideal Body Weight, Obesity Level
- **Body Composition:** Fat Mass, Muscle Mass, Body Water, Bone Mass, Subcutaneous Fat, Fat-free Body Weight, Protein Mass
- **Metabolic Indicators:** BMR (Basal Metabolic Rate), Visceral Fat, Body Age, Body Type

Measurement Results:

- After completing the measurement, basic results are displayed on the app's home screen. A detailed and accurate evaluation of individual parameters is available in the "Report" section (Figure C).

**A video tutorial with the correct procedure for using the scale is also available. The video is in English, the text description follows the currently selected app language.*

*** Heart rate measurement is not available for this model.*

Account Settings:

- In the Account -> Settings section, it is possible to adjust the basic parameters of the app. You can change the app theme (color interface), language, units of measurement, set personal goals and connect the app with Apple Health or Fitbit for data synchronization. This section also includes the option to activate a regular weighing reminder.

Troubleshooting:

1. **Bluetooth failed to connect:** Ensure the scale is on and Bluetooth is active on your phone. Check the phone's operating system version. A Bluetooth icon appears on the scale display after a successful connection. In the app, the status "Connected" appears on the main screen. If it doesn't connect automatically, go to: Account -> Device -> Add Device. Search for the Bluetooth name "WG BodyExpert" and connect manually.
2. **Body analysis is not displayed:** Measurement must be done barefoot. Ensure feet and scale surface are dry. Stand on the scale for about 10 seconds or until the value stops flashing. Leaving early may cause incomplete analysis. Check if the scale is paired with another device.
3. **Incorrect weight displayed:** Place the scale on a hard, level and stable surface. Check that the sensors (feet) on the bottom are clean. If necessary, reset the device using the reset hole located in the battery compartment. Use a thin object (e.g., a paperclip) to press and hold the button for about 6 seconds until "CLr" appears on the display.
4. **Scale won't turn on:** Check battery status and replace if necessary. Observe correct polarity.
5. **Indication symbols:**
 - Lo (Low battery): Batteries are low. Replace all 4xAA batteries.
 - Err (Overload): Scale is overloaded (exceeds 180 kg). Step off immediately.
 - C (Measurement error): Error during measurement (e.g., due to movement). Step off, wait for "0.0" and repeat.

Safety Instructions:

- Do not use if damaged or the glass is cracked. Do not attempt to disassemble or repair. All repairs should be done by an authorized service center.
- **Not for medical diagnosis.** Values are for reference only and do not replace professional medical examination.
- **Do not use if you have a pacemaker or other electronic implants.** The scale uses bioelectrical impedance, which may affect the function of implanted devices.
- Do not exceed 180 kg. Do not step on the edge or jump on the scale.
- Use only on hard, level and stable surfaces. Do not use on carpets.
- Feet and scale must be dry to prevent slipping.
- Measure at the same time of day for comparable results. Do not measure immediately after a bath or sauna.
- Pregnant women should use the weight-only function.
- Not a toy. Keep away from children. Children may use only under adult supervision.
- Use only recommended 4xAA batteries and observe polarity. Do not throw batteries into fire or short-circuit them. Dispose of batteries according to regulations.
- Clean with a slightly damp cloth; do not use aggressive chemicals.

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